

GENAXY MINI CENTRIFUGES

GENAXY Mini Centrifuge is compact and cost effective, allowing any workstation to be equipped with a "personal" mini centrifuge. It is supplied with one rotor suitable for 6x1.5/2.0ml tubes, 6x0.5ml tubes and 2x8x0.2ml PCR tubes, and two speed options (4,000 & 6,000 rpm). The electronic braking can reach quick deceleration to reduce handling time.

FEATURES

- Mini-6KS Mini-Centrifuge is compact and cost effective, allowing any workstation to be equipped with a "personal" mini centrifuge.
- It is supplied with one rotor suitable for 6x1.5/2.0ml tubes, 6x0.5ml tubes and 2x8x0.2ml PCR tubes, and two speed options (4,000 & 6,000rpm).
- The electronic braking can reach quick deceleration to reduce handling time, Composite shell, shatter-resistant.
- The electronic braking system for rapid sample processing
- Quick release rotor tub, allows for easy rotor exchange
- Two speed options: 4,000 & 6,000rpm
- Include 2 quick release interchangeable rotors, --8x1.5/2.0ml tube rotor, --4 place PCR Strip rotor



**TECHNICAL SPECIFICATION**

Model	GEM-MINI-6KS
Capacity	6x1.5/2.0ml, 6x0.5ml tubes and 2-8x0.2ml PCR Strip rotor
Two Speed options	4,000rpm & 6,000rpm
No. of UV Tubes	6 Nos.
Maximum RCF	1,200 & 2,000 x g
Dimensions	135 x 147 x 120mm (D x W x H)
Weight	0.7 Kg
Electrical	AC100~240V, 50/60Hz

AS-08031-01	6x1.5/2.0ml, 6x0.5ml tubes and 2x8x0.2ml PCR Strip rotor
--------------------	--

Genaxy Scientific Pvt. Ltd.

155-156, C-Block, Community Centre, Plot 32 to 57, Behind Janak Cinema, Janak Puri, New Delhi-110 058

Phones : +91-11-4761 9999 (100 Lines) Fax : +91-11-2553 4163

E-mail : info@genaxy.com