


# GENAXY Mini Centrifuge

## GEN-MINI-6KT

Mini-6KT Mini-Centrifuge can be supplied with different rotors for different applications. It is also offered two speed options (4,000 & 6,000 rpm). The electronic braking can reach quick deceleration to reduce handling time.



**GEN-MINI-6KT**

## Technical Specification

Model	GEN-MINI-6KT
Capacity	48 well 50ul PCR rotor or 5 x 5ml tubes rotor
Two Speed Options	4,000 rpm & 6,000 rpm
Maximum RCF	1,200 & 2,000 x g (6,000 rpm)
Dimensions (D×W×H)	135 x 147 x 120 mm
Weight	0.7 kg
Electrical	AC100~240V, 50/60Hz

## Features

Composite shell, shatter-resistant

The electronic braking system for rapid sample processing

Quick release rotor tub, allows for easy rotor exchange

Mini-6KT Two speed options: 4,000 & 6,000rpm

Patented rotor design for 48-well 50ul PCR plate Optional rotor: 4x5ml Eppendorf tube

## Accessories

S No.	Code	Description
1	AS-08101-01	48well 50ul PCR plate rotor
2	AS-08101-02	4 x 5 ml tubes rotor (Optional)
3	AS-08101-03	Adaptors for 5ml normal tubes (Optional)

### Genaxy Scientific Pvt. Ltd :

155-156, C-Block, Community Centre, Plot 32 to 57, Behind Janak Cinema, Janak Puri, New Delhi-110 058

Phones : +91-11-4761 9999 (100 Lines) Fax : +91-11-2553 4163

E-mail : info@genaxy.com