

## GENAXY Mini Plate Centrifuge

GENAXY'S Mini Plate Centrifuge is compact and economical Mini-P25 Mini-Plate Centrifuge with two-position rotor. It can quickly spin down droplets and condensation and is for use before and after thermal cycling to improve PCR yield. It accepts skirted, non-skirted and all standard PCR plates making, and it is suitable for the molecular laboratory.

### Features:-

- Quickly Spin down droplets
- Use before and after thermal cycling to increase PCR yield
- Accepts skirted, non-skirted and all standard PCR plates
- Less than 1/4 the size of most plate centrifuges



**TECHNICAL SPECIFICATION:-**

<b>Model</b>	<b>GEN-MINI-Plate</b>
Capacity	2 PCR Plates
Rotor Speed	2500 rpm
G Force	500 x g
Dimensions	20(W) x 23(D) x 19.3(H)cm
Weight	3 kg
Power	230V, 50-60HZ

**Genaxy Scientific Pvt. Ltd.**

155-156, C-Block, Community Centre, Plot 32 to 57, Behind Janak Cinema, Janak Puri, New Delhi-110 058

Phones : +91-11-4761 9999 (100 Lines) Fax : +91-11-2553 4163

E-mail : info@genaxy.com