


# MICROTITER PLATES


Genaxy Microtiter plates can be used for wide applications in the field of therapeutics, diagnostics, molecular biology, pharmaceuticals and biomedical studies.

## Unique Features

- » Manufactured with 99.9% virgin polystyrene
- » Conform to SBS/ANSI this can be used for manual and automated sample processing
- » Flat bottom and Round bottom formats with or without lid
- » High binding detachable and non-detachable plate designs are also available
- » All products are sterilized
- » Certified DNase, RNase, Human DNA and Pyrogen free



# Genaxy Microtiter Plates



Condensation rings on lid prevents cross contamination

Excellent well uniformity

Transparent round bottom surface

Transparent flat bottom surface



FOR MORE PRODUCT INFORMATION  
SCAN THE QR CODE



## Microtiter Plates

Cat. Number	Description	Qty. / Pack
GEN-MTP-96F-S	Microtitration Plates, Flat Bottom, Sterile	100 Plates
GEN-MTP-96U-S	Microtitration Plates, Round Bottom, Sterile	100 Plates
GEN-MTP-L-S	Lids for Microtiter Plates, Sterile	100 Plates
GEN-MTP-96F-L-S	Microtitration Plates with Lids, Flat Bottom, Sterile	250 Plates
GEN-MTP-96U-L-S	Microtitration Plates with Lids, Round Bottom, Sterile	250 Plates
GEN-MTP-HB-96F	96 Wells Undetachable Elisa Plate, High Binding	100 Plates
GEN-MTPS-HB-96F	96 Wells 8x12 Strips Detachable Elisa Plate, High Binding.	100 Plates

## Genaxy Scientific Pvt. Ltd.