

## Reservoir Chiller

Ice-free, non-toxic way to extend working period and preserve integrity for temperature sensitive reagents

- Reagents are kept at less than 5°C (41°F) for approximately 1 hour
- Engineered custom fit design for Heathrow Scientific 50 or 100 mL reservoirs to ensure a maximum contact for longer uniform cooling
- Eliminates the risk of contamination from messy ice bath
- Not for use with liquid nitrogen



### Symmetrical cooler design

Allows reagent reservoirs to be mounted in either direction



### Stackable

Saves storage and work space



50 mL  
Reservoir Chiller



100 mL  
Reservoir Chiller

Pack of 2

Item No.	Color	Volume mL	L x W x H in	L x W x H cm
120730	Purple	50	6.1 x 2.9 x 1.7	15.5 x 7.4 x 4.2
120731	Blue			
120732	Assorted*			
120733	Purple	100	6.2 x 9 x 1.7	15.7 x 9 x 4.2
120734	Blue			
120735	Assorted*			

\*Assorted colors (Purple and Blue)

Corporate Office :

**GENAXY SCIENTIFIC PVT. LTD.**

155-156, C-Block, Community Centre, Plot 32 to 57, Janakpuri, New Delhi - 110 058

Ph. : +91-11-4761 9999 (100Lines) Web. : [www.genaxy.com](http://www.genaxy.com) E-mail : [technical@genaxy.com](mailto:technical@genaxy.com)