



Cool Containers

Durable containers for chilling temperature-sensitive samples

Suitable for low temperature use. Ideal for ice, ice slurries and CO₂ pellets

- Made from polyurethane with metal inserts for stability
- Tapered design for easy stacking
- Close-fit lid keeps evaporation to a minimum and features a comfortable recessed knob for convenience and to save space
- Precise pour spout and large handles



Lid nests under bucket

Item No.	Color	Style	Capacity liters	L x W x H in	L x W x H cm
HS28714K	Black	■ Bucket	4	7.1 (H) x 9.2 (D)	18.1 (H) x 23.4 (D)
HS28721G	Green	■ Pan	1	10.5 x 7.7 x 4.5	26.7 x 19.6 x 11.4
HS28724B	Blue	■ Pan	4	15 x 10 x 5.5	38.1 x 25.4 x 14
HS28729R	Red	■ Pan	9	20.5 x 11.5 x 6.4	52.1 x 29.2 x 16.3

Corporate Office :

GENAXY SCIENTIFIC PVT. LTD.

155-156, C-Block, Community Centre, Plot 32 to 57, Janakpuri, New Delhi - 110 058

Ph. : +91-11-4761 9999 (100Lines) Web. : www.genaxy.com E-mail : technical@genaxy.com