

## Cryogenic Vial Canes

Organize cryogenic samples in dewars


Aluminum canes hold vials for freezer storage or immersion in LN<sub>2</sub> dewars

- Use in liquid and vapor phase; use approved vials only
- Canes bend to fit all major cryogenic vial brands



Pack of 12

Item No.	Tubes mL	No. of Vials	Length in	Length cm
HS23205	1.5/2.0	5	11.4	29
HS23206	1.5/2.0	6	11.8	30



## PVC Cryogenic Cane Sleeves

For extra security and sample safety

Clear plastic sleeves protect vials in aluminum canes during handling and storage

- Safe for cryogenic freezing
- More durable than cardboard sleeves
- Transparency allows sample identification without removing sleeve



Pack of 100

Item No.	Diameter in	Length in	Diameter cm	Length cm
HS23210	0.6	10.8	1.5	27.5



Corporate Office :

**GENAXY SCIENTIFIC PVT. LTD.**

155-156, C-Block, Community Centre, Plot 32 to 57, Janakpuri, New Delhi - 110 058  
Ph. : +91-11-4761 9999 (100Lines) Web. : [www.genaxy.com](http://www.genaxy.com) E-mail : [technical@genaxy.com](mailto:technical@genaxy.com)