

smartSlicer™

Laboratory Safety Razor Blade

Say no to dangerous sharps

- At just 0.27 mm the US Patent Pending product is as thin as a single edge industrial razor blade, but much safer
- Cleanly slice agarose gels right on the light box without scratching the surface
- Ideal for cutting soft tissue and plant matter
- Adjust the size required by simply cutting with scissors
- Made from food grade PET



Pack of 100

Item No.	L x W x H in	
120259	0.8 x 1.5 x 0.01	19 x 38 x 0.3
120260	1 x 2 x 0.01	25 x 50 x 0.3

smartSlicer™ is a trademark of LevGo Inc.

Corporate Office :

GENAXY SCIENTIFIC PVT. LTD.

155-156, C-Block, Community Centre, Plot 32 to 57, Janakpuri, New Delhi - 110 058

Ph. : +91-11-4761 9999 (100Lines) Web. : www.genaxy.com E-mail :

technical@genaxy.com