



**Weighted Stand**  
for bench top use

**Lock Counter Off**  
to allow use of the marker  
without registering more counts

**Hand Tally Adapter**  
allows you to use the eCount  
as a digital tally counter

## eCount Colony Counter

Combines a multi-function electronic counter with a Sharpie® Pen for marking to prevent missing or double-counting colonies

- Light, ergonomic body is balanced for easy manipulation
- Select between an audible beep or an LED flash to acknowledge count
- Store up to 32 separate counts and display the total of all counts
- Count up/down verification
- Quick and easy to change between different colored markers to check different colony counts in the same view



Battery and Pen not included in warranty

US Patent No. D564,385  
EP Registered Design No. 705751-0001

Item No.	L x W in	L x W cm
120000	7.4 x 1.2*	17 x 3*
120228	CR1225 Battery	

\*Measurements shown are with Sharpie® Pen inserted

### Package Includes

eCount Colony Counter	Hand Tally Adapter
Black Sharpie® Pen	Weighted Stand
Long-life Button Battery	



Corporate Office :

**GENAXY SCIENTIFIC PVT. LTD.**

155-156, C-Block, Community Centre, Plot 32 to 57, Janakpuri, New Delhi - 110 058  
Ph. : +91-11-4761 9999 (100Lines) Web. : www.genaxy.com E-mail : technical@genaxy.com