

Inoculating Needle/Loops

Combination provides choice to spread or pick colonies and cultures

Transfer bacterial cultures without risk of contamination

- Use loop end to uniformly streak a gel surface
- Use needle end for removal of colonies
- Choose between flexible polypropylene or the more rigid ABS versions
- Lot-to-lot traceability

STERILE EO



Item No.	Color	Material	Loop μ L	L x D in	L x D cm	Packaging 1,000 per Box	Sterile
HS81121A	Blue ■	Flexible PP	1	7.9 x .2	20 x .4	50 x 20pk	✓
HS81121B		Rigid ABS					
HS81121C	Blue ■	Flexible PP	10	7.9 x .2	20 x .4	50 x 20pk	✓
HS81121D		Rigid ABS					

Corporate Office :

GENAXY SCIENTIFIC PVT. LTD.

155-156, C-Block, Community Centre, Plot 32 to 57, Janakpuri, New Delhi - 110 058

Ph. : +91-11-4761 9999 (100Lines) Web. : www.genaxy.com E-mail : technical@genaxy.com