

Hold standard slides (76 x 26 mm and 3 x 1 in)

Cardboard slide folders offer a convenient solution for safely working, transporting and short-term storing of slides

Organize & Protect

- Index on the front cover corresponds with the number grid inside for identification
- Recessed slide compartments and inside covers for specimen protection
- Slides cannot be accessed when cover is closed reducing contamination

Sturdy Construction

- Compartment holder is made from a double layer of densely pressed cardboard
- Cover is single layer cardboard
- Front protection flaps are attached with durable tape

Item No.	Color	Capacity	L x W x H in	L x W x H cm
HS9902	Black	■		
HS9913	Blue	■		
HS9914	Green	■	8.2 x 13.4 x 0.4	20.8 x 34.1 x 0.9
120196	Red	■		
120197	Yellow	■		
HS9930	Black	■	8.2 x 12.9 x 0.4	28.1 x 32.8 x 1



Corporate Office :

GENAXY SCIENTIFIC PVT. LTD.

155-156, C-Block, Community Centre, Plot 32 to 57, Janakpuri, New Delhi - 110 058
Ph. : +91-11-4761 9999 (100Lines) Web. : www.genaxy.com E-mail : technical@genaxy.com