



Slide Holders

Hold slides during evaluation, cleaning and drying

- Made from chemical-resistant polyethylene making the unit resist stains and easy to clean
- Solid, strong construction with rubber feet for bench-top stability
- Angled slide holder slots allow easy access to individual slides and viewing of slide labels



Item No.	Color	Capacity	L x W x H <i>in</i>	L x W x H <i>cm</i>
HS15981A	Natural <input type="checkbox"/>	40-Place	11.7 x 4.1 x 0.9	29.8 x 10.5 x 2.4
HS15981B	Natural <input type="checkbox"/>	78-Place	11.8 x 7 x 0.8	29.9 x 17.8 x 2.1
HS15981C	Natural <input type="checkbox"/>	15-Place	6 x 4 x 0.8	15.3 x 10.2 x 2.1