

Slide Staining Tray Set

Reduce possible cross-contamination and messy work areas

This rugged staining tray is an economical alternative to metal versions

- Made from a polypropylene and polyethylene blend for stability
- Stain, rinse, and dry slides on a single working tray
- A clean and fresh working surface area every time
- Up to eight slides fit comfortably onto the base and the deep well holds up to 38 mL
- Compact size and recessed handles for easy transport from work area to sink
- Pack includes 4 trays and 1 lid



Item No.	Color	Capacity	L x W x H <i>in</i>	L x W x H <i>cm</i>
HS15951A	Black ■	8-Place	12 x 5 x 1.2	30.5 x 12.7 x 3

Corporate Office :

GENAXY SCIENTIFIC PVT. LTD.

155-156, C-Block, Community Centre, Plot 32 to 57, Janakpuri, New Delhi - 110 058

Ph. : +91-11-4761 9999 (100Lines) Web. : www.genaxy.com E-mail : technical@genaxy.com