

## Slide Tray Clasps

Effortlessly locate, retrieve and keep slide folders organized

- Assorted color clasps are ideal for color coding identification and to keep slide folder flaps securely closed
- Clasps can be labeled with a permanent marker or a label can be attached

### Slide Tray Clasps

Pack of 100

Item No.	Color	L x W x H in	L x W x H cm
120653	Assorted* <input type="checkbox"/> <input type="checkbox"/> <input type="checkbox"/> <input type="checkbox"/> <input type="checkbox"/>	2.2 x 1.7 x 0.6	5.6 x 4.4 x 1.6
120654	Clear <input type="checkbox"/>		

\*Assorted colors (Clear, Grey, Blue, Green, Red)

