

Cell Spreaders

In microbiology, a cell spreader or plate spreader is a tool used to smoothly spread cells and bacteria on a culture plate, such as a petri dish.



Delta™

- Made from ABS
- Flexible rippled handle and domed spreading bar to help reduce agar tears

STERILE R



L-Shape

- Made from ABS
- Lot-to-lot traceability
- 10 Pack comes in a zip bag

STERILE R



Stainless Steel

- Flame-sterilizable
- Reusable



Item No.	Shape	L x W x H in	L x W x H cm	Packaging	Sterile
HS8151	Delta™	6.8 x 1.5 x .2	17.2 x 3.7 x .4	500/Single Wrap	✓
HS8171A	L-Shape	5.7 x 1.5 x .2	14.6 x 3.8 x .4	500/Single Wrap	✓
HS8171B			500, 50 x 10 Pack	✓	
HS86655	Stainless Steel	8 x 1.3 x .2	20.3 x 3.3 x .4	1/Single Wrap	