



Sliding Lock
provides secure storage with easily identified lock and unlock symbols.

Slide Rail Design

ensures a secure fit allowing easy placement and removal of slides

Centrally Positioned Numbering

enables simple slide organization with easy viewing of location numbers



100-Place Premium Plus Slide Box

Modern design for maximum functionality

Optimal Organization

- Organize by “project” and “date” molded in on lid
- Inventory card index corresponds to slide numbers and is removable for trouble-free writing on flat surface
- Vibrant translucent colors are easy to identify and enable coding options

Functional Design

- Flexible storage, store vertically or horizontally to efficiently use space
- Molded feet allow secure horizontal stacking of multiple units
- Unique side storage can accommodate a pen, desiccant or other small items needed for working
- Hinged lid smoothly opens 180° sitting flat on surfaces and prevents misplacement

Safe & Secure

- EVA foam lines the bottom of the box helping to hold slides securely in place and protect against damage
- Stylish wave surface texture is easy to hold for safe sample transport

100-Place Premium Plus Slide Box

Pack of 5

Item No.	Color	Item No.	Color
120576	Blue	120579	Red
120577	Orange	120592	Purple
120578	Green	120593	Assorted*

*Assorted colors (Purple, Blue, Green, Red, Orange)

Lining	Lock	L x W x H in	L x W x H cm
EVA Foam	Sliding Latch	8.2 x 7.7 x 1.3	20.7 x 19.5 x 3.2

Replacement Inventory Card

Pack of 10

Item No.	L x H in	L x H cm
120618	8.5 x 11	21.6 x 28

Corporate Office :

GENAXY SCIENTIFIC PVT. LTD.

155-156, C-Block, Community Centre, Plot 32 to 57, Janakpuri, New Delhi - 110 058

Ph. : +91-11-4761 9999 (100Lines) Web. : www.genaxy.com E-mail : technical@genaxy.com