

## Standard Microscope Slide Boxes

Constructed of durable polypropylene



Slide boxes hold either 25 or 100 standard slides (76 x 26 mm and 3 x 1 in)

### Improved Identification

- Molded slot numbers correspond directly to the slide inventory card
- Organize with "Project" and "Date" pre-molded labeling area

### Secure & Stackable

- Hinged lids are held securely closed with a latch lock
- Molded feet secure stacking for easy transport and storage

### 25-Place Slide Box

Item No.	Color	Lining	Lock	L x W x H in (cm)
HS15990A	Blue	None	Latch Lock	5.6 x 3.6 x 1.2 (14.1 x 9.2 x 3)
HS15990B	Red			
HS15990C	White			
120740	25-Place Replacement Inventory Card, Pack of 10			5.3 x 3.1 (13.4 x 8)

### 100-Place Slide Box

Item No.	Color	Lining	Lock	L x W x H in (cm)
HS15991A	Blue	None	Latch Lock	8.3 x 6.7 x 1.5 (21 x 16.9 x 3.7)
HS15991B	Red			
HS15991C	White			
120741	100-Place Replacement Inventory Card, Pack of 10			3.9 x 6.1 (10 x 15.5)

## Low Temperature Microscope Slide Boxes

Store in ultra-low temperatures with confidence



Holds standard slides (76 x 26 mm and 3 x 1 in)

### Durable Construction

- Made of polycarbonate to store in ultra-low temperatures
- Uniform grooved slide slots separate slides preventing slide contact
- Molded feet enable secure stacking for easy transport and storage

### Durable Construction

- Low temperature resistant numbered foam lining corresponds with identification sheet
- Organize with "Project" and "Date" pre-molded labeling area



### Low Temperature Microscope Slide Boxes

Item No.	Color	Capacity	Lock	L x W x H in (cm)
120316	Blue	25-Place	Friction Fit	5.6 x 3.5 x 1.4 (14.1 x 8.8 x 3.5)
HS15988H	Purple	100-Place	Stainless Steel Clasp	8.2 x 6.9 x 1.3 (20.8 x 17.5 x 3.4)