



AEROSOL FILTER BARRIER GENTIPS

OPTIMAL PROTECTION HOLDS BACK
— AEROSOL & CONTAMINATION —

Utilizing a unique manufacturing process, Genaxy Scientific has developed an innovative series of pipette tips that feature ultra smooth surfaces. Genaxy Filter Barrier Tips blocks aerosol by coalescing radioactive, chemicals, or biological contaminants before they can come in contact with pipettor shaft. This effectively eliminates unwanted signals induced by residual RNA/DNA transferred in trace amounts by aerosol cross contamination. They also act as a safety barrier against accidental over pipetting of fluids without risk of sample lock-up.

Unique Features

- » Made of 99.9% virgin Polypropylene conforming to US FDA 21 CFR.
- » Manufactured under ISO Class 8 clean room facility
- » Filters are made up of Polyethylene matrix with optimized airflow that is 100% inert, hydrophobic and pure. This protects against aerosol bypass and sample carryover, eliminates the threat of leachables, extractables and delivers the maximum sample recovery
- » Internal ring design ensures maximum compatibility with popular pipettes.
- » Diamond polished injection mould gives the smoothest internal surface.
- » Zero Fluid Retention (ZFR) technology offers maximum fluid recovery.
- » Accurate graduation marks for precise liquid aspiration.
- » Perfect fine point orifice for the highest levels of accuracy and reproducibility.
- » Tips are DNase/RNase, Human DNA, Heavy Metal, Pyrogen, Phthalate and PCR inhibitor free.
- »



Aerosol Filter Barrier GenTips

Maximize Accuracy & Precision



Genfit ensure the maximum compatibility with popular pipettes, so reducing the number of the tip styles & brand on your lab bench

Genaxy filter barrier tips block aerosols by coalescing radioactive, chemical, or biological contaminants before they can come in contact with pipette shaft. this effectively eliminates unwanted signals induced by residual RNA/DNA transferred in trace amounts by aerosol cross-contamination. they also act as a safety barrier against accidental over-pipetting of fluids without risk of sample lock-up

Injection mold core pins are diamond polished to ensure the smoothest internal surface, so reducing liquid hand-up as well as providing the highest transparency to improve content visibility

Zero fluid retention technology offers unparalleled maximum fluid recovery to save precious sample & reagent. the food colouring additive in the tip shown adjacent, visually demonstrates the superior low retention capabilities of national tips over other manufacturers

Perfect fine point orifice for the highest level of accuracy & reproducibility



FOR MORE PRODUCT INFORMATION
SCAN THE QR CODE



Aerosol/Filter Barrier Pipette Tips

Cat. Number	Description	Qty. / Pack
GEN-FB-10	0.5-10µl Aerosol Barrier Gentip, Natural, Bulk, Non-Sterile	1000 tips
GEN-FB-10-S	0.5-10µl Aerosol Barrier Gentip, Natural, Bulk, Sterile	1000 tips
GEN-FB-10-L	0.5-10µl Aerosol Barrier Gentip, Natural, Bulk, Low Retention, Non-Sterile	1000 tips
GEN-FB-10-L-S	0.5-10µl Aerosol Barrier Gentip, Natural, Bulk, Low Retention, Sterile	1000 tips
GEN-FB-10-R-S	0.5-10µl Aerosol Barrier Gentip, Natural, Racked, Sterile	10 racks
GEN-FB-10-L-R-S	0.5-10µl Aerosol Barrier Gentip, Natural, Racked, Low Retention, Sterile	10 racks
GEN-FB-10XL	0.5-20µl Aerosol Barrier Gentip, Natural, Bulk, Non-Sterile	1000 tips
GEN-FB-10XL-S	0.5-20µl Aerosol Barrier Gentip, Natural, Bulk, Sterile	1000 tips
GEN-FB-10XL-L	0.5-20µl Aerosol Barrier Gentip, Natural, Bulk, Low Retention, Non-Sterile	1000 tips
GEN-FB-10XL-L-S	0.5-20µl Aerosol Barrier Gentip, Natural, Bulk, Low Retention, Sterile	1000 tips
GEN-FB-10XL-R-S	0.5-20µl Aerosol Barrier Gentip, Natural, Racked, Sterile	10 racks
GEN-FB-10XL-L-R-S	0.5-20µl Aerosol Barrier Gentip, Natural, Racked, Low Retention, Sterile	10 racks
GEN-FB-20	0.5-20µl Aerosol Barrier Gentip, Natural, Bulk, Non-Sterile	1000 tips
GEN-FB-20-S	0.5-20µl Aerosol Barrier Gentip, Natural, Bulk, Sterile	1000 tips
GEN-FB-20-L	0.5-20µl Aerosol Barrier Gentip, Natural, Bulk, Low Retention, Non-Sterile	1000 tips
GEN-FB-20-L-S	0.5-20µl Aerosol Barrier Gentip, Natural, Bulk, Low Retention, Sterile	1000 tips
GEN-FB-20-R-S	0.5-20µl Aerosol Barrier Gentip, Natural, Racked, Sterile	10 racks
GEN-FB-20-L-R-S	0.5-20µl Aerosol Barrier Gentip, Natural, Racked, Low Retention, Sterile	10 racks
GEN-FB-50	0.5-50µl Aerosol Barrier Gentip, Natural, Bulk, Non-Sterile	1000 tips
GEN-FB-50-S	0.5-50µl Aerosol Barrier Gentip, Natural, Bulk, Sterile	1000 tips
GEN-FB-50-L	0.5-50µl Aerosol Barrier Gentip, Natural, Bulk, Low Retention, Non-Sterile	1000 tips
GEN-FB-50-L-S	0.5-50µl Aerosol Barrier Gentip, Natural, Bulk, Low Retention, Sterile	1000 tips
GEN-FB-50-R-S	0.5-50µl Aerosol Barrier Gentip, Natural, Racked, Sterile	10 racks
GEN-FB-50-L-R-S	0.5-50µl Aerosol Barrier Gentip, Natural, Racked, Low Retention, Sterile	10 racks
GEN-FB-100	1-100µl Aerosol Barrier Gentip, Natural, Bulk, Non-Sterile	1000 tips
GEN-FB-100-S	1-100µl Aerosol Barrier Gentip, Natural, Bulk, Sterile	1000 tips
GEN-FB-100-L	1-100µl Aerosol Barrier Gentip, Natural, Bulk, Low Retention, Non-Sterile	1000 tips
GEN-FB-100-L-S	1-100µl Aerosol Barrier Gentip, Natural, Bulk, Low Retention, Sterile	1000 tips
GEN-FB-100-R-S	1-100µl Aerosol Barrier Gentip, Natural, Racked, Sterile	10 racks
GEN-FB-100-L-R-S	1-100µl Aerosol Barrier Gentip, Natural, Racked, Low Retention, Sterile	10 racks
GEN-FB-200	1-200µl Aerosol Barrier Gentip, Natural, Bulk, Non-Sterile	1000 tips

Aerosol/Filter Barrier Pipette Tips

Cat. Number	Description	Qty. / Pack
GEN-FB-200-S	1-200µl Aerosol Barrier Gentip, Natural, Bulk, Sterile	1000 tips
GEN-FB-200-L	1-200µl Aerosol Barrier Gentip, Natural, Bulk, Low Retention, Non-Sterile	1000 tips
GEN-FB-200-L-S	1-200µl Aerosol Barrier Gentip, Natural, Bulk, Low Retention, Sterile	1000 tips
GEN-FB-200-R-S	1-200µl Aerosol Barrier Gentip, Natural, Racked, Sterile	10 racks
GEN-FB-200-L-R-S	1-200µl Aerosol Barrier Gentip, Natural, Racked, Low Retention, Sterile	10 racks
GEN-FB-300	1-300µl Aerosol Barrier Gentip, Natural, Bulk, Non-Sterile	1000 tips
GEN-FB-300-L	1-300µl Aerosol Barrier Gentip, Natural, Bulk, Low Retention, Non-Sterile.	1000 tips
GEN-FB-300-R-S	1-300µl Aerosol Barrier Gentip, Natural, Racked, Sterile.	10 racks
GEN-FB-300-L-R-S	1-300µl Aerosol Barrier Gentip, Natural, Racked, Low Retention, Sterile.	10 racks
GEN-FB-1000	100-1000µl Aerosol Barrier Gentip, Natural, Bulk, Non-Sterile.	1000 tips
GEN-FB-1000-S	100-1000µl Aerosol Barrier Gentip, Natural, Bulk, Sterile.	1000 tips
GEN-FB-1000-L	100-1000µl Aerosol Barrier Gentip, Natural, Bulk, Low Retention, Non-Sterile.	1000 tips
GEN-FB-1000-L-S	100-1000µl Aerosol Barrier Gentip, Natural, Bulk, Low Retention, Sterile.	1000 tips
GEN-FB-1000-R-S	100-1000µl Aerosol Barrier Gentip, Natural, Racked, Sterile.	6 racks
GEN-FB-1000-L-R-S	100-1000µl Aerosol Barrier Gentip, Natural, Racked, Low Retention, Sterile.	6 racks
GEN-FBG-1000	100-1000µl Aerosol Barrier Gentip, Natural, Bulk, Graduated, Non-Sterile.	1000 tips
GEN-FBG-1000-S	100-1000µl Aerosol Barrier Gentip, Natural, Bulk, Graduated, Sterile.	1000 tips
GEN-FBG-1000-L	100-1000µl Aerosol Barrier Gentip, Natural, Bulk, Graduated, Low Retention, Non-Sterile.	1000 tips
GEN-FBG-1000-L-S	100-1000µl Aerosol Barrier Gentip, Natural, Bulk, Graduated, Low Retention, Sterile.	1000 tips
GEN-FBG-1000-R-S	100-1000µl Aerosol Barrier Gentip, Natural, Racked, Graduated, Sterile.	10 Racks
GEN-FBG-1000-L-R-S	100-1000µl Aerosol Barrier Gentip, Natural, Racked, Graduated, Low Retention, Sterile	10 Racks
GEN-FBG-1250	100-1250µl Aerosol Barrier Gentip, Natural, Bulk, Graduated, Non-Sterile.	1000 tips
GEN-FBG-1250-S	100-1250µl Aerosol Barrier Gentip, Natural, Bulk, Graduated, Sterile.	1000 tips
GEN-FBG-1250-L	100-1250µl Aerosol Barrier Gentip, Natural, Bulk, Graduated, Low Retention, Non-Sterile.	1000 tips
GEN-FBG-1250-L-S	100-1250µl Aerosol Barrier Gentip, Natural, Bulk, Graduated, Low Retention, Sterile.	1000 tips

Aerosol/Filter Barrier Pipette Tips

Cat. Number	Description	Qty. / Pack
GEN-FBG-1250-R-S	100-1250µl Aerosol Barrier Gentip, Natural, Racked, Graduated, Sterile.	8 Racks
GEN-FBG-1250-L-R-S	100-1250µl Aerosol Barrier Gentip, Natural, Racked, Graduated, Low Retention, Sterile.	8 Racks
GEN-FB-1250XL	100-1250µl Aerosol Barrier Gentip, Natural, Bulk, Non-Sterile.	1000 tips
GEN-FB-1250XL-S	100-1250µl Aerosol Barrier Gentip, Natural, Bulk, Sterile.	1000 tips
GEN-FB-1250XL-L	100-1250µl Aerosol Barrier Gentip, Natural, Bulk, Low Retention, Non-Sterile.	1000 tips
GEN-FB-1250XL-L-S	100-1250µl Aerosol Barrier Gentip, Natural, Bulk, Low Retention, Sterile.	1000 tips
GEN-FB-1250XL-R-S	100-1250µl Aerosol Barrier Gentip, Natural, Racked, Sterile.	6 Racks
GEN-FB-1250XL-L-R-S	100-1250µl Aerosol Barrier Gentip, Natural, Racked, Low Retention, Sterile.	6 Racks
GEN-DFB-10-S	0.5-10µl Dual Aerosol Barrier Gentip, Natural, Bulk, Sterile.	1000 tips
GEN-DFB-10XL	0.5-10µl Dual Aerosol Barrier Gentip, Natural, Bulk, Non-Sterile.	1000 tips
GEN-DFB-10XL-S	0.5-10µl Dual Aerosol Barrier Gentip, Natural, Bulk, Sterile.	1000 tips
GEN-DFB-10-R-S	0.5-10µl Dual Aerosol Barrier Gentip, Natural, Racked, Sterile.	10 racks
GEN-DFB-10XL-R-S	0.5-10µl Dual Aerosol Barrier Gentip, Natural, Racked, Sterile.	10 racks
GEN-DFB-200	1-200µl Dual Aerosol Barrier Gentip, Natural, Bulk, Non-Sterile.	1000 tips
GEN-DFB-200-S	1-200µl Dual Aerosol Barrier Gentip, Natural, Bulk, Sterile.	1000 tips
GEN-DFB-200-R-S	1-200µl Dual Aerosol Barrier Gentip, Natural, Racked, Sterile.	10 racks
GEN-DFB-1000	100-1000µl Dual Aerosol Barrier Gentip, Natural, Bulk, Non-Sterile.	1000 tips
GEN-DFB-1000-S	100-1000µl Dual Aerosol Barrier Gentip, Natural, Bulk, Sterile.	1000 tips
GEN-DFB-1000-R-S	100-1000µl Dual Aerosol Barrier Gentip, Natural, Racked, Sterile.	6 racks

Genaxy Scientific Pvt. Ltd.