

TISSUE & CELL DNA/RNA Kits

FavorPrep™ Tissue and Cell DNA Kits

For Extraction of Genomic DNA from animal tissues, fixed tissues, bacteria (G+ and G-) yeast, dried blood or cultured cells without using organic solvents.

- ❑ The procedure takes about 20 minutes after lysis. Yields high quality DNA up to 50µg from mini kit.
- ❑ The eluted DNA is ready to use in different applications like PCR, AFLP, RFLP, Southern Blotting and Real Time PCR.
- ❑ These kits are time saving, versatile and safe.
- ❑ The kits are available in Mini & 96 well format



CATALOG NO.	Description	Kit Size
FATGK 001	FavorPrep™ Tissue Genomic DNA Extraction Mini Kit	50 Preps
FATGK 001-1	FavorPrep™ Tissue Genomic DNA Extraction Mini Kit	100 Preps
FATGK 001-2	FavorPrep™ Tissue Genomic DNA Extraction Mini Kit	300 Preps
FATGM 001	Tissue Genomic DNA Extraction MicroElute Kit with Proteinase K Powder	50 Preps
FATGM 001-1	Tissue Genomic DNA Extraction MicroElute Kit with Proteinase K Powder	100 Preps
High throughput 96 Well Format		
FADWE 96004	FavorPrep™ 96 Well Genomic DNA Extraction Kit	96 x 4 Plates

FavorPrep™ Tissue and Cell RNA Kits

For Extraction of total RNA from animal tissues, cultured cells bacteria and yeast without using organic solvents.

- ❑ The procedure takes about 60 minutes and Yields high quality RNA up to 40µg from mini kit, 1.5mg from Maxi Kit and about 60µg from cultured cell.
- ❑ These kits are time saving and safe.
- ❑ The RNA is highly pure with OD 260/280: 1.9-2.1
- ❑ The kits are available in Mini, Maxi& 96 well format.

Besides total RNA kit for miRNA Isolation is also available. This Kit is designed for purification of small RNA with minimal contamination from large RNA molecules or genomic DNA. The purified miRNA can be used for gene expression, cancer markers and developmental control. Expected yield of miRNA is about 2µg.



Catalog No.	Description	Kit Size
FATRK 001	FavorPrep™ Tissue total RNA Purification Mini Kit	50 Preps
EATRK 001-1	FavorPrep™ Tissue total RNA Purification Mini Kit	100 Preps
FATRK 001-2	FavorPrep™ Tissue total RNA Purification Mini Kit	300 Preps
FATRK 003	FavorPrep™ Tissue total RNA Purification Maxi Kit	10 Preps
FATRK 003-1	FavorPrep™ Tissue total RNA Purification Maxi Kit	24 Preps
FAMIK 001	FavorPrep™ miRNA Isolation Kit	100 Preps
FAMIK 002	FavorPrep™ miRNA Isolation Kit	50 Preps
FATRE 96004	FavorPrep™ 96 Well Total RNA Kit	96 x 4 Plates

