

NITRILE GLOVES

- 5Mil (5.5g) powder free Nitrile Examination gloves.
- Exceptional puncture and chemical resistance with excellent barrier quality.
- Free from any natural rubber protein.
- Ambidextrous with textured fingers and beaded cuff.
- ISO 9001 certified and FDA registered.
- ASTM D6319, EN 455 PART 1,2,3 & 4 compliant.
- Ice blue (bluple) in color.
- Registered with CDSCO (Central Drugs Standard Control Organization)

Cat. Number	Description	Qty. / Pack
GEN-NXG-S	Nitrile Exam Gloves, Bluple, Small Size	50 pair
GEN-NXG-M	Nitrile Exam Gloves, Bluple, Medium Size	50 pair
GEN-NXG-L	Nitrile Exam Gloves, Bluple, Large Size	50 pair

Genaxy Scientific Pvt. Ltd.