



SCREW CAP MICRO TUBES

Screw cap tubes are ideal for long term use in a variety of situations. The cap and O-ring are designed to be the perfect complement to the screw tube thus ensuring a perfect seal.

- + Tube/Cap = virgin polypropylene
- + O-ring = silicon rubber (Please loosen the cap and autoclave)
- + Centrifugal force: 20,000RCF (physiologic(al) saline 1.8mL/75minutes), Flat top (recordable)



0.5ml Screw Cap Tubes

Catalog	Description	Style	Cap Color	Type	Qty./Pack
GEN-ST-050-C	Screw Cap Tubes with Caps	Conical	Natural	Non-Sterilized	500 X 2
GEN-ST-050-C-S	Screw Cap Tubes with Caps	Conical	Natural	Pre-Sterilized	500 X 2
GEN-ST-050-A	Screw Cap Tubes with Caps	Conical	Assorted*	Non-Sterilized	500 X 2
GEN-ST-050-A-S	Screw Cap Tubes with Caps	Conical	Assorted*	Pre-Sterilized	500 X 2
GEN-ST-050-COLOR*	Screw Cap Tubes with Caps	Conical	W/P/R/G/Y/B	Non-Sterilized	500 X 2
GEN-ST-050-X	Screw Cap Amber Tubes with Caps	Conical	Brown	Non-Sterilized	500 X 2
GEN-ST-050-X-S	Screw Cap Amber Tubes with Caps	Conical	Brown	Pre-Sterilized	500 X 2
GEN-ST-050-SS-C	Screw Cap Tubes with Caps	Self Standing	Natural	Non-Sterilized	500 X 2
GEN-ST-050-SS-C-S	Screw Cap Tubes with Caps	Self Standing	Natural	Pre-Sterilized	500 X 2
GEN-ST-050-SS-A	Screw Cap Tubes with Caps	Self Standing	Assorted	Non-Sterilized	500 X 2
GEN-ST-050-SS-A-S	Screw Cap Tubes with Caps	Self Standing	Assorted	Pre-Sterilized	500 X 2
GEN-ST-050-SS-COLOR*	Screw Cap Tubes with Caps	Self Standing	W/P/R/G/Y/B	Non-Sterilized	500 X 2
GEN-ST-050-SS-X	Screw Cap Amber Tubes with Caps	Self Standing	Brown	Non-Sterilized	500 X 2
GEN-ST-050-SS-X-S	Screw Cap Amber Tubes with Caps	Self Standing	Brown	Pre-Sterilized	500 X 2



*For the Color of your choice, Simply Filling the "COLOR" with desired color code. (W-White, P-Pink, B-Blue, G-Green, O-Orange, R-Red, Y-Yellow)

1.5ml Screw Cap Tubes

Catalog	Description	Style	Cap Color	Type	Qty./Pack
GEN-ST-150-C	Screw Cap Tubes with Caps	Conical	Natural	Non-Sterilized	500 X 2
GEN-ST-150-C-S	Screw Cap Tubes with Caps	Conical	Natural	Pre-Sterilized	500 X 2
GEN-ST-150-A	Screw Cap Tubes with Caps	Conical	Assorted*	Non-Sterilized	500 X 2
GEN-ST-150-A-S	Screw Cap Tubes with Caps	Conical	Assorted*	Pre-Sterilized	500 X 2
GEN-ST-150-COLOR*	Screw Cap Tubes with Caps	Conical	W/P/R/G/Y/B	Non-Sterilized	500 X 2
GEN-ST-150-X	Screw Cap Amber Tubes with Caps	Conical	Brown	Non-Sterilized	500 X 2
GEN-ST-150-X-S	Screw Cap Amber Tubes with Caps	Conical	Brown	Pre-Sterilized	500 X 2
GEN-ST-150-SS-C	Screw Cap Tubes with Caps	Self Standing	Natural	Non-Sterilized	500 X 2
GEN-ST-150-SS-C-S	Screw Cap Tubes with Caps	Self Standing	Natural	Pre-Sterilized	500 X 2
GEN-ST-150-SS-A	Screw Cap Tubes with Caps	Self Standing	Assorted	Non-Sterilized	500 X 2
GEN-ST-150-SS-A-S	Screw Cap Tubes with Caps	Self Standing	Assorted	Pre-Sterilized	500 X 2
GEN-ST-150-SS-COLOR*	Screw Cap Tubes with Caps	Self Standing	W/P/R/G/Y/B	Non-Sterilized	500 X 2
GEN-ST-150-SS-X	Screw Cap Amber Tubes with Caps	Self Standing	Brown	Non-Sterilized	500 X 2
GEN-ST-150-SS-X-S	Screw Cap Amber Tubes with Caps	Self Standing	Brown	Pre-Sterilized	500 X 2



*For the Color of your choice, Simply Filling the "COLOR" with desired color code. (W-White, P-Pink, B-Blue, G-Green, O-Orange, R-Red, Y-Yellow)



2.0ml Screw Cap Tubes

Catalog	Description	Style	Cap Color	Type	Qty./Pack
GEN-ST-200-C	Screw Cap Tubes with Caps	Conical	Natural	Non-Sterilized	500 X 2
GEN-ST-200-C-S	Screw Cap Tubes with Caps	Conical	Natural	Pre-Sterilized	500 X 2
GEN-ST-200-A	Screw Cap Tubes with Caps	Conical	Assorted*	Non-Sterilized	500 X 2
GEN-ST-200-A-S	Screw Cap Tubes with Caps	Conical	Assorted*	Pre-Sterilized	500 X 2
GEN-ST-200-COLOR*	Screw Cap Tubes with Caps	Conical	W/P/R/G/Y/B	Non-Sterilized	500 X 2
GEN-ST-200-X	Screw Cap Amber Tubes with Caps	Conical	Brown	Non-Sterilized	500 X 2
GEN-ST-200-X-S	Screw Cap Amber Tubes with Caps	Conical	Brown	Pre-Sterilized	500 X 2
GEN-ST-200-SS-C	Screw Cap Tubes with Caps	Self Standing	Natural	Non-Sterilized	500 X 2
GEN-ST-200-SS-C-S	Screw Cap Tubes with Caps	Self Standing	Natural	Pre-Sterilized	500 X 2
GEN-ST-200-SS-A	Screw Cap Tubes with Caps	Self Standing	Assorted	Non-Sterilized	500 X 2
GEN-ST-200-SS-A-S	Screw Cap Tubes with Caps	Self Standing	Assorted	Pre-Sterilized	500 X 2
GEN-ST-200-SS-COLOR*	Screw Cap Tubes with Caps	Self Standing	W/P/R/G/Y/B	Non-Sterilized	500 X 2
GEN-ST-200-SS-X	Screw Cap Amber Tubes with Caps	Self Standing	Brown	Non-Sterilized	500 X 2
GEN-ST-200-SS-X-S	Screw Cap Amber Tubes with Caps	Self Standing	Brown	Pre-Sterilized	500 X 2

*For the Color of your choice, Simply Filling the "COLOR" with desired color code. (W-White, P-Pink, B-Blue, G-Green, O-Orange, R-Red, Y-Yellow)

TRANSPORT TUBES

Our durable and reliable Transport Tube and Cap provide a 'leak-proof' and reliable seal for your critical sample storage. Both Tube and Cap are manufactured of virgin polypropylene giving our product the ability to be autoclaved.

- + Graduation marks at 1.0 mL intervals
- + Large, frosted writing area for sample identification
- + Secure fit between tube and cap provides a leak resistant seal



Catalog	Description	Volume	Tube Color	Type	Qty./Pack
GEN-ST-5ml	Transport Tubes	5ml	Clear	Non-Sterilized	500
GEN-ST-5ml-S	Transport Tubes	5ml	Clear	ETO Sterile	500
GEN-ST-12ml	Transport Tubes	12ml	Clear	Non-Sterilized	500
GEN-ST-12ml-S	Transport Tubes	12ml	Clear	ETO Sterile	500

RIA VIAL

These Polypropylene/Polystyrene tubes are generally used in RIA, Coagulation and Bacteriology. The Polypropylene tubes are clear and can withstand 3000 RCF during Centrifugation while the transparent Polystyrene tubes can withstand 1400 RCF. The PP tubes can accept common acids, alkalies & solvents at room temperature while the PS tubes can accept mild bases, weak acids but not organic solvents.



Catalog	Description	Ø x H (mm)	Material	Volume	Qty./Pack
GEN-1275-PS	Ria Vial Tubes	12 X 75	PS	5ml	500
GEN-1275-PP	Ria Vial Tubes	12 X 75	PP	5ml	500
GEN-1275-PP-X	Ria Vial Amber Tubes	12 X 75	PP	5ml	500
GEN-12-CAP	Ria Vial Caps	12	HDPE	-	500
GEN-12-CAP-X	Ria Vial Amber Caps	12	HDPE	-	500