



100 WELL STORAGE RACK

Made from high grade polypropylene Genaxy storage racks are useful for safe and secure storage of microtube and cryogenic vials. Easy to open hinged lid with snap clasp. Alphanumeric coding helps to identify and locate the samples. Vibrant colors are available, special black rack for light sensitive samples. All racks are autoclavable. Dimension: 14.1 x 15.1 x 5.7 cm ; Well dia. 1.23 cm

Cat. No.	Product Description	Pack Size
GEN-CSR-100	100-Well PP cryobox, Natural	5
GEN-CSR-100-A	100-Well PP cryobox, Assorted Racks	5

Genaxy Scientific Pvt. Ltd .