

# CARBOYS WITH TUBULATION

## Carboys With Tubulation Polypropylene



- Option for Medical Grade PP
- Medical Grade conforming to USP Class VI, US FDA CFR.
- Meets Chemical Registration of TSCA (USA) & DSL (Canada).
- Autoclavable

Cat. Number	Description	Volume	Type/Style	Qty. / Pack
GEN-CBTB-PP-10L	CARBOY	10 Ltr.	TUBULATION, PP	1
GEN-CBTB-PP-20L	CARBOY	20 Ltr.	TUBULATION, PP	1
GEN-CBTB-PP-10L-MG	CARBOY	10 Ltr.	TUBULATION, PP Medical Grade	1
GEN-CBTB-PP-20L-MG	CARBOY	20 Ltr.	TUBULATION, PP Medical Grade	1

## Carboys with Tubulation LDPE



- Low Density Polyethylene (LDPE).
- Option for Medical Grade LDPE
- Medical Grade conforming to USP Class VI, US FDA CFR.
- Meets Chemical Registration of TSCA (USA) & DSL (Canada).

Cat. Number	Description	Volume	Type/Style	Qty. / Pack
GEN-CBTB-LDPE-10L	CARBOY	10 Ltr.	TUBULATION, LDPE	1
GEN-CBTB-LDPE-20L	CARBOY	20 Ltr.	TUBULATION, LDPE	1
GEN-CBTB-LDPE-10L-MG	CARBOY	10 Ltr.	TUBULATION, LDPE Medical Grade	1
GEN-CBTB-LDPE-20L-MG	CARBOY	20 Ltr.	TUBULATION, LDPE Medical Grade	1

**Genaxy Scientific Pvt. Ltd.**