


# PCR PLATES

NOVASBIO multi-well PCR Plates, with their ultra-thin wall design, ensure rapid and efficient heat transfer between tube contents and cycler block.

## 96 well, No Skirt

- Broad cycler compatibility
- Available in Low Profile 0.1ml version for FAST blocks
- Single notch corner at A12




**Cycler: Compatibility** Analytik/Jena/Bioметра qTOWER<sup>3</sup> / G / touch, Topical, T1 Thermocycler, Tgradient, Tone, Tadvanced, Trobot 96, SpeedCycler2, Bio-rad MyiQ<sup>™</sup>, iCycler<sup>™</sup>, iCycler<sup>™</sup> IQ / IQ 4 / IQ 5, MyCycler<sup>™</sup>, T100, C1000 Touch, S1000, Eppendorf MasterCycler<sup>®</sup> ep realplex, ep / epgradient / Pro / Pro S / Cycler:nexus / nexus gradient / nexus SX1

CATALOG NO.	VOLUME	DESCRIPTION	PACK	CASE
PCR96-FLT-C	0.2ml	PCR 96 Well Plates, No-Skirt, Clear	10 Plates	5 Packs

## 96 well, Low profile, No Skirt

- Broad cycler compatibility
- Low Profile wells with reduced volume (0.1ml)
- Available in white for increased autofluorescence
- Dual notch corner at A12, H12




**Cycler: Compatibility** Agilent AriaMx, Analytik/Jena/Bioметра Trobot 96, SpeedCycler2 (SP, SPR), Bio-rad CFX96 Touch, CFX96 Touch Deep Well, CFX connect, C1000 Touch, S1000, Eppendorf MasterCycler<sup>®</sup> ep realplex, ep/epgradient / Pro / Pro S / Cycler:nexus / nexus gradient / nexus SX1

CATALOG NO.	VOLUME	DESCRIPTION	PACK	CASE
PCR96LP-FLT-C	0.2ml	PCR 96 Well Plates, No-Skirt, Clear	10 Plates	10 Packs
PCR96LP-FLT-W	0.2ml	PCR 96 Well Plates, No-Skirt, White	10 Plates	10 Packs

## 96 well, Half Skirt

- Half skirt style provides increased rigidity
- Traditional profile, flat-top plate
- Skirt allows for labeling and barcodes
- High contrast alphanumeric lettering for easy well identification
- Single notch corner at A12



**Cycler: Compatibility** Agilent Mx3000P<sup>™</sup>, Mx3005P<sup>™</sup>, Mx4000P<sup>™</sup>, Analytik/Jena/Bioметра qTOWER<sup>3</sup>/G/touch, Topical, Flexcycler2, T1 Thermocycler, Tgradient, Tone, Tadvanced, TProfessional Gradient/XL, Trobot 96, SpeedCycler2 Bio-rad MyiQ<sup>™</sup>, iCycler<sup>™</sup> IQ / IQ 4 / IQ 5, MyCycler<sup>™</sup>, T100, C1000 Touch, S1000 Eppendorf Master Cycler<sup>®</sup> ep realplex, ep/epgradient / Pro / Pro S / Cycler:nexus / nexus gradient / nexus SX1

CATALOG NO.	VOLUME	DESCRIPTION	PACK	CASE
PCR96-HS-C	0.2ml	PCR 96 Well Plates, Clear, Half Skirt	10 Plates	5 Packs




Manufactured from virgin, medical grade polypropylene and certified free from RNase, DNase, Human DNA and pyrogens, these plates offer exceptional compatibility across the major cycler brands.


## 96 well, Half Skirt, ABI compatible

- Elevated rim around plate
- Fits ABI Standard Cyclers and Sequencers
- Standard volume wells (0.2ml)
- Single notch corner at A12

For SeptaMat,  
see page 10



Cycler: **ABI QuantStudio™ 3 / 5 / 6 / 7 / 12K, ViiA7™**  
Compatibility: 7000, 7300, 7500, 7700, 7900 HT **GeneAmp® 2700 / 2720 / 9600 / 9700, Veriti, Proflex, Simplicamp 3100, 3130XL, 3500, 3500XL, 3730, 3730XL**



CATALOG NO.	VOLUME	DESCRIPTION	PACK	CASE
PCR96-AB-C	0.2ml	PCR 96 Well Plates, Clear, Half Skirt	10 Plates	5 Packs

## 96 well, Half Skirt, Low profile, ABI compatible

- Elevated rim around plate
- Fits ABI FAST Cyclers
- Low Profile wells with reduced volume (0.1ml)
- Single notch corner at A1

Cycler: **Applied Biosystems QuantStudio™ 3/5/6/7/12K, ViiA7™**  
Compatibility: 7500 FAST, 7900 HT FAST, StepOne™, StepOne™ Plus, GeneAmp® 9800 FAST, Veriti FAST




CATALOG NO.	VOLUME	DESCRIPTION	PACK	CASE
PCR96-LP-AB-C	0.1ml	PCR 96 Well Plates, Clear, Half Skirt	10 Plates	10 Packs

## 384 well, Full Skirt

- Automation-friendly
- Broad compatibility with 384-block thermal cyclers
- Clear version with single notch at A24
- White version designed for Roche® Lightcycler® 480 with dual notch corner at A24, P24

Cycler: **ABI, GeneAmp® 9700, Veriti, Proflex, Multiblock system 3100, 3130XL, 3500, 3500XL, 3730, 3730XL, AnalytikJena/Biomtra Flexcyler2, T1Thermocycler, Advanced, TProfessional, Trobot 96, Bio-rad CFX384 Touch, C1000 Touch, S1000**



CATALOG NO.	VOLUME	DESCRIPTION	PACK	CASE
PCR384-C	50µl	PCR 384 Well Plates, Clear, Full Skirt	10 Plates	5 Packs
PCR384-W	50µl	PCR 384 Well Plates, White, Full Skirt	10 Plates	5 Packs