

Hard Shell PCR Plates



Roche® compatible
HP96-LPSSR-CC

Rigid polycarbonate frame eliminates warping during PCR

Ultra clear thin-walled polypropylene wells for optimum real-time PCR and regular PCR

Recommended for low volume PCR reactions

Compatible with Bio-Rad®CFX Touch systems, Agilent AriaMx, Roche® LightCycler®, Analytik Jena qTOWER²/ qTOWER Touch and Techne Quantica

CATALOG NO.	VOLUME	DESCRIPTION	PACK	CASE
HP96-LPFS-BKC	0.1ml	Hard Shell 96-well PCR Plates, Low Profile, Skirted, Black frame, Clear wells	10 Plates	5 Packs
HP96-LPFS-WW	0.1ml	Hard Shell 96-well PCR Plates, Low Profile, Skirted, White frame, White wells	10 Plates	5 Packs
HP96-LPSSR-CC	0.1ml	Hard shell 96-well PCR Plates, Semi-skirted. Low Profile. Clear frame. Clear wells	10 Plates	5 Packs
HP96-LPSSR-WW	0.1ml	Hard shell 96-well PCR Plates, Semi-skirted. Low Profile. White frame. White wells	10 Plates	5 Packs