

# SAMPLE CONTAINERS

## SAMPLE CONTAINERS



- Ultraclear , polypropylene sample containers
- Thick walls for sturdiness
- LOCK SYSTEM: The external thread of the container has a small triple bridge type structure which while closing the cap comes in contact with the teethered structure present in the internal thread of the cap. This results in tight and secure locking of the cap.
- Clear embossed graduation marks for measuring
- Sterile containers are ETO sterile and come in individually wrapped packaging

Cat. Number	Description	Volume	Packaging	Type/Style	Qty. / Pack
GEN-SC-30	Sample Container	30ml	Bulk	Non-sterile	500
GEN-SC-30-S	Sample Container	30ml	Individual Pack	ETO Sterile	300
GEN-SC-50	Sample Container	50ml	Bulk	Non-sterile	200
GEN-SC-50-S	Sample Container	50ml	Individual Pack	ETO Sterile	100
GEN-SC-70	Sample Container	70ml	Bulk	Non-sterile	200
GEN-SC-70-S	Sample Container	70ml	Individual Pack	ETO Sterile	100
GEN-SC-100	Sample Container	100ml	Bulk	Non-sterile	200
GEN-SC-100-S	Sample Container	100ml	Individual Pack	ETO Sterile	100

**Genaxy Scientific Pvt. Ltd.**