

# GENAXY<sup>®</sup>

Customer Satisfaction Through Business Ethics



Brushless DC motor



Inbuilt counter balance  
(No walking issues)



Microprocessor based  
closed loop control



## G-MIX MINI VORTEX

## G-MIX

Mini Vortex Mixer G-MIX is a compact unit equipped with a brushless DC motor that provide fast and efficient mixing with minimal vibration. Delivers a speed of 3000 Rpm, and orbital diameter of 4mm, and has a pulse mode function to increase the vortexing efficiency.

- ◀ Adjustable speed from 1000 to 3000 RPM & timer setting
- ◀ Analog display to shows speeds
- ◀ Unique pulse mode function – User enabled
- ◀ (On/Off Switch) with Power On LED
- ◀ Orbital diameter of 4 mm
- ◀ Can be used in hoods & cold rooms
- ◀ Various attachments are available for different capacities of tubes.



## Technical Specification

Catalogue No.	G-MIX
Max Speed	3000 RPM
Orbital diameter	4 mm
Motor Type	Brushless DC
Digital Display	Yes
Over all Dimensions (W X H)	130 x 170 mm
Speed Setting	1000-3000 RPM
Timer Setting	NA
Pulse mode	Yes
Max Load Capacity	500 gm
Weight	1.5 kg
Attachments	Round carrier / microtube & cup as Standard Attachments

## Ordering Information

Catalogue No.	G-MIX
G-MIX	Analog Vortex Mixer – Max RPM of 3000

