



UNIVERSAL FIT GENTIPS


LOW RETENTION

HIGH PRECISION

GENTip universal tip comes with a bevelled or fine point orifice. Bevelled orifice which gives shorter contact area while aspirating and dispensing, results in No Last Drop Issue and No False Results. While the standard tips, without bevelled orifice, comes with larger surface area, which reduces the control, results in last drop issues and no accuracy. GenTips are available in bulk packs, racked form and individually packed. Bulk tips are bagged, and heat sealed in "puncture resistant" laminated nylon-base dispenser packs with super wide gusset at bottom to stand neatly upright. Stack rack options are also available in 10ul, 200ul and 1000ul. All racks are sterile and individually wrapped.

Unique Features

- » Made of 99.9% medical grade virgin Polypropylene which meets USP class IV requirements. High purity and chemical resistance.
- » PP material is free of any additives, plasticizers, or colorants to prevent sample contamination.
- » Conforms to USFDA 21 CFR and manufactured under ISO Class 8 clean room facility.
- » ISO 9001:2015, 13485:2016 & 14001:2015 certified.
- » Unique unigrip design with 3 internal rings creating a vacuum allowed onto the pipette shaft eliminating sideways dispensing
- » Use of high-quality diamond polished injection mould for smooth internal surface and reduced liquid hang-up.
- » Zero Fluid Retention (ZFR) technology offers maximum fluid recovery and minimize liquid adhesion, ensuring complete sample dispensing.
- » Accurate graduation marks provide visual confirmation that the correct liquid volume is being aspirated or dispensed.
- » Tips are free from particles and fibers, ensuring high transparency and better visualization of samples.
- » Touch centrum design gives universal fit to all popular brands like Eppendorf, Finnpipette, Pipetman (Gilson) etc.
- » Tips are DNase/RNase, Human DNA, Heavy metal, Pyrogen(Endotoxin), Phthalate and PCR inhibitor free (Reports are available on request).
- » Autoclavable for sterilization at 121°C for 15 minutes.
Sterile tips are gamma sterile.




Genaxy Universal Fit GenTips



Beveled Oriŷces For Total Sample Recovery

Blade Prevents Droplets & Increases Precision



Molded In Exclusive 99.9% Virgin Polypropylene Resin Which Results In Exceptional Sample Delivery & Improved Assay Performance

High-quality diamond polished injection mould for smooth internal surface and reduced liquid hang-up



Better Sealability With Soft But Touch Centrum Design



FOR MORE PRODUCT INFORMATION
SCAN THE QR CODE



Bulk & Racked Pipette Tips

Cat. Number	Description	Qty. / Pack
GEN-UT-10-C	0.5-10µl Universal Fit Gentip, Natural, Bulk, Non-Sterile.	1000 tips
GEN-UT-10-L-C	0.5-10µl Universal Fit Gentip, Natural, Bulk, Low Retention, Non-Sterile.	1000 tips
GEN-UT-10-C-R	0.5-10µl Universal Fit Gentip Natural, Racked, Non-Sterile.	10 racks
GEN-UT-10-C-R-S	0.5-10µl Universal Fit Gentip, Natural, Racked, Sterile.	10 racks
GEN-UT-10-C-STK	0.5-10µl Universal Fit Gentip, Natural, Reload Stacks, Non-Sterile.	960 Tips
GEN-UT-10-L-C-STK	0.5-10µl Universal Fit Gentip, Natural, Reload Stacks, Low Retention, Non-Sterile.	10 racks
GEN-FB-10XL	0.5-20µl Aerosol Barrier Gentip, Natural, Bulk, Non-Sterile	1000 tips
GEN-FB-10XL-S	0.5-20µl Aerosol Barrier Gentip, Natural, Bulk, Sterile	1000 tips
GEN-FB-10XL-L	0.5-20µl Aerosol Barrier Gentip, Natural, Bulk, Low Retention, Non-Sterile	1000 tips
GEN-FB-10XL-L-S	0.5-20µl Aerosol Barrier Gentip, Natural, Bulk, Low Retention, Sterile	1000 tips
GEN-FB-10XL-R-S	0.5-20µl Aerosol Barrier Gentip, Natural, Racked, Sterile	10 racks
GEN-FB-10XL-L-R-S	0.5-20µl Aerosol Barrier Gentip, Natural, Racked, Low Retention, Sterile	10 racks
GEN-FB-20	0.5-20µl Aerosol Barrier Gentip, Natural, Bulk, Non-Sterile	1000 tips
GEN-FB-20-S	0.5-20µl Aerosol Barrier Gentip, Natural, Bulk, Sterile	1000 tips
GEN-FB-20-L	0.5-20µl Aerosol Barrier Gentip, Natural, Bulk, Low Retention, Non-Sterile	1000 tips
GEN-FB-20-L-S	0.5-20µl Aerosol Barrier Gentip, Natural, Bulk, Low Retention, Sterile	1000 tips
GEN-FB-20-R-S	0.5-20µl Aerosol Barrier Gentip, Natural, Racked, Sterile	10 racks
GEN-FB-20-L-R-S	0.5-20µl Aerosol Barrier Gentip, Natural, Racked, Low Retention, Sterile	10 racks
GEN-FB-50	0.5-50µl Aerosol Barrier Gentip, Natural, Bulk, Non-Sterile	1000 tips
GEN-FB-50-S	0.5-50µl Aerosol Barrier Gentip, Natural, Bulk, Sterile	1000 tips
GEN-FB-50-L	0.5-50µl Aerosol Barrier Gentip, Natural, Bulk, Low Retention, Non-Sterile	1000 tips
GEN-FB-50-L-S	0.5-50µl Aerosol Barrier Gentip, Natural, Bulk, Low Retention, Sterile	1000 tips
GEN-FB-50-R-S	0.5-50µl Aerosol Barrier Gentip, Natural, Racked, Sterile	10 racks
GEN-FB-50-L-R-S	0.5-50µl Aerosol Barrier Gentip, Natural, Racked, Low Retention, Sterile	10 racks
GEN-FB-100	1-100µl Aerosol Barrier Gentip, Natural, Bulk, Non-Sterile	1000 tips
GEN-FB-100-S	1-100µl Aerosol Barrier Gentip, Natural, Bulk, Sterile	1000 tips
GEN-FB-100-L	1-100µl Aerosol Barrier Gentip, Natural, Bulk, Low Retention, Non-Sterile	1000 tips
GEN-FB-100-L-S	1-100µl Aerosol Barrier Gentip, Natural, Bulk, Low Retention, Sterile	1000 tips
GEN-FB-100-R-S	1-100µl Aerosol Barrier Gentip, Natural, Racked, Sterile	10 racks
GEN-FB-100-L-R-S	1-100µl Aerosol Barrier Gentip, Natural, Racked, Low Retention, Sterile	10 racks
GEN-FB-200	1-200µl Aerosol Barrier Gentip, Natural, Bulk, Non-Sterile	1000 tips
GEN-FB-200-S	1-200µl Aerosol Barrier Gentip, Natural, Bulk, Sterile	1000 tips
GEN-FB-200-L	1-200µl Aerosol Barrier Gentip, Natural, Bulk, Low Retention, Non-Sterile	1000 tips
GEN-FB-200-L-S	1-200µl Aerosol Barrier Gentip, Natural, Bulk, Low Retention, Sterile	1000 tips
GEN-FB-200-R-S	1-200µl Aerosol Barrier Gentip, Natural, Racked, Sterile	10 racks
GEN-FB-200-L-R-S	1-200µl Aerosol Barrier Gentip, Natural, Racked, Low Retention, Sterile	10 racks
GEN-FB-300	1-300µl Aerosol Barrier Gentip, Natural, Bulk, Non-Sterile	1000 tips
GEN-FB-300-L	1-300µl Aerosol Barrier Gentip, Natural, Bulk, Low Retention, Non-Sterile.	1000 tips

Bulk & Racked Pipette Tips

Cat. Number	Description	Qty. / Pack
GEN-UT-200-C-L-R-S	1-200µl Universal Fit Gentip, Natural, Racked, Bevelled, Low Retention, Sterile	10 racks
GEN-UTG-200-C-R	1-200µl Universal Fit Gentip, Graduated, Natural, Racked, Bevelled, Non-Sterile.	10 racks
GEN-UTG-200-Y-R	1-200µl Universal Fit Gentip, Graduated, Yellow, Racked, Bevelled, Non-Sterile.	10 racks
GEN-UTG-200-C-R-S	1-200µl Universal Fit Gentip, Graduated, Natural, Racked, Bevelled, Sterile.	10 racks
GEN-UTG-200-Y-R-S	1-200µl Universal Fit Gentip, Graduated, Yellow, Racked, Bevelled, Sterile.	10 racks
GEN-UTG-200-C-L-R	1-200µl Universal Fit Gentip, Graduated, Natural, Racked, Bevelled, Low Retention, Non-Sterile	10 racks
GEN-UTG-200-C-L-R-S	1-200µl Universal Fit Gentip, Graduated, Natural, Racked, Bevelled, Low Retention, Sterile	10 racks
GEN-UT-200-C-STK	1-200µl Universal Fit Gentip, Natural, Reload Stacks, Bevelled, Non-Sterile	960 tips
GEN-UT-200-Y-STK	1-200µl Universal Fit Gentip, Yellow, Reload Stacks, Bevelled, Non-Sterile	960 tips
GEN-UT-200-LC-STK	1-200µl Universal Fit Gentip, Natural, Reload Stacks, Bevelled, Low Retention, Non-Sterile	10 racks
GEN-UT-300-C	1-300µl Universal Fit Gentip, Natural, Bulk, Non-Sterile	1000 tips
GEN-TT-300-C	1-300µl Universal Fit Gentip, Natural, Bulk, Non-Sterile, Finn timer Compatible	1000 tips
GEN-UT-300-C-R	1-300µl Universal Fit Gentip, Natural, Racked Non-Sterile	10 racks
GEN-UT-300-C-R-S	1-300µl Universal Fit Gentip, Natural Racked, Sterile	10 racks
GEN-TT-300-C-R	1-300µl Universal Fit Gentip, Natural, Racked, Non-Sterile, Finn timer Compatible	10 racks
GEN-TT-300-C-R-S	1-300µl Universal Fit Gentip, Natural, Racked, Sterile, Finn timer Compatible	10 racks
GEN-UT-1000-C	100-1000µl Universal Fit Gentip, Natural, Bulk, Bevelled, Non-Sterile	500 tips
GEN-UT-1000-B	100-1000µl Universal Fit Gentip, Blue, Bulk, Bevelled, Non-Sterile	500 tips
GEN-UTG-1000-C	100-1200µl Universal Fit Gentip, Natural, Bulk, Graduated, Non-Sterile (83mm)	500 tips
GEN-UTG-1000-B	100-1200µl Universal Fit Gentip, Blue, Bulk, Graduated, Non-Sterile (83mm)	500 tips
GEN-UTG-1250-C	100-1250µl Universal Fit Gentip, Natural, Bulk, Graduated, Non-Sterile (104mm)	500 tips
GEN-TE-1000-B	100-1000µl Universal Fit Gentip, Blue, Bulk, Eppendorf Style, Non-Sterile	10x500 tips
GEN-UT-1000-L-C	100-1000µl Universal Fit Gentip, Natural, Bulk, Bevelled, Low Retention, Non-Sterile	500 tips

Bulk & Racked Pipette Tips

Cat. Number	Description	Qty. / Pack
GEN-UTG-1000-L-C	100-1200µl Universal Fit Gentip, Natural, Bulk, Bevelled, Low Retention, Non-Sterile	500 tips 10 racks
GEN-UT-1000-C-R	100-1000µl Universal Fit Gentip, Natural, Racked, Bevelled, Non-Sterile	10 racks
GEN-UT-1000-B-R	100-1000µl Universal Fit Gentip, Blue, Racked, Bevelled, Non-Sterile	
GEN-UTG-1000-C-R	100-1200µl Universal Fit Gentip, Natural, Racked, Graduated, Non-Sterile (83mm)	10 racks
GEN-UTG-1250-C-R	100-1250µl Universal Fit Gentip, Natural, Racked, Graduated, Non-Sterile (104mm)	8 racks
GEN-UTG-1000-C-STK	1-1000µl Universal Fit Gentip, Natural, Reload Stacks, Graduated, Non-Sterile	10 racks
GEN-UTG-1000-LC-STK	1-1000µl Universal Fit Gentip, Natural, Reload Stacks, Graduated, Low Retention, Non-Sterile	10 racks
GEN-UT-1000-C-R-S	100-1000µl Universal Fit Gentip, Natural, Racked, Bevelled, Sterile	10 racks
GEN-UT-1000-B-R-S	100-1000µl Universal Fit Gentip, Blue, Racked, Bevelled, Sterile	10 racks
GEN-UTG-1000-C-R-S	100-1200µl Universal Fit Gentip, Natural, Racked, Graduated, Sterile (83mm)	10 racks
GEN-UTG-1250-C-R-S	100-1250µl Universal Fit Gentip, Natural, Racked, Graduated, Sterilized (104mm)	8 racks
GEN-UT-1000-C-L-R	100-1000µl Universal Fit Gentip, Natural, Racked, Bevelled, Low Retention, Non-Sterile	10 racks
GEN-UT-1000-C-L-R-S	100-1000µl Universal Fit Gentip, Natural, Racked, Bevelled, Low Retention, Sterile	10 racks
GEN-UT-1250XL-C	100-1200µl Universal Fit Gentip, Natural, Bulk, Non-Sterile (120mm), Compatible with 1ml Pipettors	500 tips
GEN-UT-1250XL-C-R	100-1200µl Universal Fit Gentip, Natural, Racked, Non-Sterile.(120mm), Compatible with 1ml Pipettors	6 racks
GEN-UT-1250XL-C-R-S	100-1200µl Universal Fit Gentip, Natural, Racked, Sterile.(120mm), Compatible with 1ml Pipettors	6 racks
GEN-UT-5000-C	1-5ml Universal Fit Gentip, Natural, Bulk, Gilson Style, Non-Sterile	250 tips
GEN-ET-5000-C	1-5ml Universal Fit Gentip, Natural, Bulk, Eppendorf Style, Non-Sterile	300 tips
GEN-TT-5000-C	1-5ml Universal Fit Gentip, Natural, Bulk, Thermo Finntip Style, Non-Sterile	300 tips
GEN-UT-10ML-C	1-10ml Universal Fit Gentip, Natural, Bulk, Gilson Style, Non-Sterile	200 tips

Individually Wrapped Pipette Tips

Cat. Number	Description	Qty. / Pack
GEN-Y-200-IW	200µl Universal Fit Gentip, Yellow, Individual Peel O° , Plastic Wrap, Sterile	500 tips
GEN-C-200-IW	200µl Universal Fit Gentip, Clear, Individual Peel O° , Plastic Wrap, Sterile	500 tips
GEN-C-200-IW-PP	200µl Universal Fit Gentip, Clear, Individual Peel O° , Paper Plastic Wrap, Sterile	5 x 400 tips
GEN-C-300-IW	10µl Universal Fit Gentip, Clear, Individual Peel O° , Plastic Wrap, Sterile.	500 tips
GEN-B-1000-IW	1000µl Universal Fit Gentip, Blue, Individual Peel O° , Plastic Wrap, Sterile.	500 tips
GEN-C-1000-IW	1000µl Universal Fit Gentip, Clear, Individual Peel O° , Plastic Wrap, Sterile.	500 tips
GEN-C-1000-IW-PP	1000µl Universal Fit Gentip, Clear, Individual Peel O° , Paper Plastic Wrap, Sterile.	5 x 200 tips

Empty Racks for Pipette Tips

GEN-UT-10-ER	Empty Racks for GenTip 10µl	10
GEN-UT-200-ER	Empty Racks for GenTip 200µl	10
GEN-UT-300-ER	Empty Racks for GenTip 300µl	10
GEN-ET-200-ER	Empty Racks for Eppendorf Style Tips,200µl	10
GEN-UT-1000-ER	Empty Racks for GenTip 1000µl	10
GEN-UTG-1000-ER	Empty Racks for GenTip 1200µl	10
GEN-UTG-1250-ER	Empty Racks for GenTip 1250µl	8
GEN-UTG-1250XL-ER	Empty Racks for GenTip 1250µl	6
GEN-ET-1000-ER	Empty Racks for Eppendorf Style Tips, 1000µl	10
GEN-UT-0510-ER	Empty Racks for 5ml & 10ml Tips.(Accepts Gilson Style Tips)	Each

Genaxy Scientific Pvt. Ltd.

155-156, C-Block, Community Centre, Plot 32 to 57, Behind Janak Cinema, Janak Puri, New Delhi-110 058
 Phones : +91-11-4761 9999 (100 Lines) Fax : +91-11-2553 4163, E-mail : info@genaxy.com