


“just keep spinning”



CENTRIFUGE TUBES

ELITE PERFORMANCE

— AT —

EXCEPTIONAL VALUE




GENAXY Brand Centrifuge Tubes are designed to meet the needs Of Researchers. This clear polypropylene Centrifuge tubes comes with the advanced High Density Polyethylene caps packed in high quality ziplock bags. These tubes are available in 15 ml, 25ml, 30 ml and 50 ml volume with amber (15ml, 50ml Tubes Only) and clear color options. Genaxy also provide these tubes in sterile, non-sterile, racked and individually wrapped formats.

Unique Features

- » The tubes are manufactured from medical grade polypropylene which meets USP class IV requirements.
- » Conforms to US FDA 21 CFR and manufactured Under ISO Class 8 Clean Room Facility.
- » ISO 9001:2015, 13485:2016 & 14001:2015 certified.
- » High spin rates with maximum RCF ratings of up to 18,000 x g (15ml & 50ml Tubes).
- » 25ml tubes will spin up to 12,500xg RCF.
- » Temperature ranges from -80°C to 121°C. Autoclavable.
- » Easy Grip HDPE caps have an integral seal ring design for secure closures that eliminate seepage when used as recommended.
- » Clear, Black, accurate printed graduations with unique Silk Printing technology resisting solvents like Acetone, Methanol etc.
- » Cap can be closed with quick 2-3 turns allows single hand operation for closing and opening of the tube.
- » Large white writing area for sample labeling.
- » Ultra-smooth conical bottoms, which allow for 100% visibility and no sharp "gates" to cause cuts and scratches.
- » Sterile tubes are radiation sterilized which ensures the non-cytotoxicity.
- » Certified DNase/RNase, Heavy-metal, Human DNA, Phthalate and pyrogen free
- » Sterile tubes are Gamma irradiated to SAL 10-6.



Genaxy Centrifuge Tubes



Flat Cap For Labeling


Easy Grip Cap with integral seal ring for easy gripping and one hand operation



Writing Area For Sample Identification



Extremely Accurate Printed Graduations With unique silk printing



Ultra Smooth Bottom For 100% Visibility- No Sharp 'Gates' To Cause Cuts And Scratches



FOR MORE PRODUCT INFORMATION
SCAN THE QR CODE





Cat. Number	Description	Qty. / Pack
GEN-CT-15ML	15ml Centrifuge Tube with Cap, Conical, Clear, Non-Sterile	25 X 20
GEN-CT-15ML-S	15ml Centrifuge Tube with Cap, Conical, Clear, Sterile	25 X 20
GEN-CT-15ML-R	15ml Centrifuge Tube with Cap, Conical, Racked, Clear, Non-Sterile	25 X 20
GEN-CT-15ML-R-S	15ml Centrifuge Tube with Cap, Conical, Racked, Clear, Sterile	25 X 20
GEN-CT-15ML-X	15ml Centrifuge Tube with Cap, Conical, Amber, Non-Sterile	25 X 20
GEN-CT-15ML-X-S	15ml Centrifuge Tube with Cap Conical, Amber, Sterile	25 X 20
GEN-CT-15ML-X-R	15ml Centrifuge Tube with Cap, Conical, Racked, Clear, Non-Sterile	25 X 20
GEN-CT-15ML-X-R-S	15ml Centrifuge Tube with Cap, Conical, Racked, Clear, Sterile	25 X 20
GEN-CT-25ML-S	25ml Centrifuge Tube with Cap, Conical, Clear, Sterile	25 X 8
GEN-CT-25ML-R-S	25ml Centrifuge Tube with Cap, Conical, Racked, Clear, Non-Sterile	25 X 8
GEN-CT-50ML	50ml Centrifuge Tube with Cap, Conical, Clear, Non-Sterile	25 X 20
GEN-CT-50ML-S	50ml Centrifuge Tube with Cap, Conical, Clear, Pre- Sterilized	25 X 20
GEN-CT-50ML-R	50ml Centrifuge Tube with Cap, Conical, Racked, Clear, Non-Sterile	12 X 25
GEN-CT-50ML-R-S	50ml Centrifuge Tube with Cap, Conical, Racked, Clear, Sterile	12 X 25
GEN-CT-50ML-X	50ml Centrifuge Tube with Cap, Conical, Amber, Non-Sterile	25 X 20
GEN-CT-50ML-X-S	50ml Centrifuge Tube with Cap, Conical, Amber, Sterile	25 X 20
GEN-CT-50ML-X-R	50ml Centrifuge Tube with Cap, Conical, Racked, Amber, Non-Sterile.	12 X 25
GEN-CT-50ML-X-R-S	50ml Centrifuge Tube with Cap, Conical, Racked, Amber, Sterile.	12 X 25
GEN-CT-50ML-SS	50ml Centrifuge Tube with Cap, Self Standing, Clear, Non-Sterile	25 X 20
GEN-CT-50ML-SS-S	50ml Centrifuge Tube with Cap, Self Standing, Clear, Sterile	25 X 20
GEN-CT-30ML-SS-S	30ml Centrifuge Tube with Cap, Self Standing, Clear, Sterile	500
GEN-CT-15ML-S-IW	15ml Centrifuge Tube with Cap, Conical, Clear, Individually Wrapped, Sterile	250
GEN-CT-50ML-S-IW	50ml Centrifuge Tube with Cap, Conical, Clear, Individually Wrapped, Sterile	250
GEN-CT-10-UP	Centrifuge Tube with Plug Cap, 10ml	500
GEN-CT-10-US	Centrifuge Tube with Screw Cap, 10ml	500
GEN-CT-50-UP	Centrifuge Tube with Plug Cap, 50ml	400
GEN-CT-50-US	Centrifuge Tube with Screw Cap, 50ml	400

Genaxy Scientific Pvt. Ltd.

155-156, C-Block, Community Centre, Plot 32 to 57, Behind Janak Cinema, Janak Puri, New Delhi-110 058
 Phones : +91-11-4761 9999 (100 Lines) Fax : +91-11-2553 4163, E-mail : info@genaxy.com