

Droplet™ Sample Storage Tray

Sample organization for the busy lab

Hard wearing storage trays that can be used for specimen containers with the insert provided, or for general lab products without the insert

- Lightweight and sturdy to hold up to 24 full specimen containers (Specimen containers not included)
- Stackable with other trays to support as many as four high
- Can be easily cleaned with common disinfectants and germicides
- Bin and tray are made from polystyrene



Use for storage or transport of laboratory products



Removable Inserts
allow outer shell to be used for storage or transport of laboratory products


Integrated Label Holder

Item No.	Color	Wells	Well Dia. mm	L x W x H in	L x W x H cm
120172	Blue/White	24	58	18.9 x 13 x 3.3	48 x 33.1 x 8.4

Coated Wire Rack - Urine Containers

Durable racks for transport or sample handling

- Chemically-resistant, HDPE-coated steel wire rack
- Fits bottles with 55-58 mm diameters
- Mesh bottom support



Item No.	Color	Tube Array	Wells	Well Dia. mm	L x W x H in	L x W x H cm
120091	Blue	6 x 4	24	58	15.5 x 10.3 x 2	39.3 x 26.2 x 5