



# TRANSPORT TUBES

*MOVE SAFELY*


— *WITH* —

*US*

Genaxy durable and reliable Transport Tube and Cap provide a 'leak-proof' and reliable seal for your critical sample storage. Both Tube and Cap are manufactured of virgin polypropylene giving our product the ability to be autoclaved.

## Unique Features

- » Manufactured with High grade 99.9% Virgin polypropylene under ISO class 8 clean room facility
- » Conforms to US FDA 21 CFR
- » Leak-proof and reliable seal for critical sample storage
- » Graduation marks at 1 ml intervals
- » Large, frosted writing area for sample labelling
- » Sterile or Non-sterile forms are available
- » Sterile Tubes are either ETO sterile or Gamma sterile (both options are given)
- » 5ml and 12ml volumes are available
- » Certified DNase/RNase, Human DNA and Phthalate free



# Genaxy Transport Tubes



Leak-proof and reliable seal ensures safe transportation of samples



Accurate graduation for visual confirmation



Frosted area for writing



Skirted bottom for better stability on benchtop



FOR MORE PRODUCT INFORMATION  
SCAN THE QR CODE



## Transport Tubes

Cat. Number	Description	Qty. / Pack
GEN-ST-5ml	5ml Transport Tubes, Clear, Non-Sterile	100 X 5
GEN-ST-5ml-S	5ml Transport Tubes, Clear, Sterile (ETO Sterile)	100 X 5
GEN-ST-5ml-GS	5ml Transport Tubes, Clear, Sterile (GAMMA Sterile)	25 X 20
GEN-ST-12ml	12ml Transport Tubes, Clear, Non-Sterile	100 X 5
GEN-ST-12ml-S	12ml Transport Tubes, Clear, Sterile	100 X 5
GEN-ST-12ml-GS	12ml Transport Tubes, Clear, Sterile (GAMMA Sterile)	25 X 20

**Genaxy Scientific Pvt. Ltd.**

155-156, C-Block, Community Centre, Plot 32 to 57, Behind Janak Cinema, Janak Puri, New Delhi-110 058  
Phones : +91-11-4761 9999 (100 Lines) Fax : +91-11-2553 4163, E-mail : [info@genaxy.com](mailto:info@genaxy.com)