

GENAXY ROUND BOTTOM TEST TUBE (FACS)



- 5ml and 14 mL polystyrene round bottom test tubes have a 1,400 RCF rating.
- The tubes come with polyethylene dual-position snap caps, which offer both vented and fully closed positions.
- RNase/DNase-free.
- Widely used and referenced in laboratory protocols
- Dual position snap-cap prevents sample loss with ease of handling
- Sterile, Non-pyrogenicity tested to less than 0.1 EU/ml

Cat. Number	Description	Volume	Size (ØxH) mm	Type/Style	Qty. / Pack
GEN-14100-PS	Test Tubes, Round Bottom	14ml	17 X 100	PS	500
GEN-14100-PS-IW	Test Tubes, Round Bottom, Individually Wrapped	14ml	17 X 100	PS	500
GEN-1275-SCPS	Test Tubes, Round Bottom	5ml	12 X 75	PS	1000
GEN-1275-SCPS-IW	Test Tubes, Round Bottom, Individually Wrapped	5ml	12 X 75	PS	500

Genaxy Scientific Pvt. Ltd.

155-156, C-Block, Community Centre, Plot 32 to 57, Behind Janak Cinema, Janak Puri, New Delhi-110 058
 Phones : +91-11-4761 9999 (100 Lines) Fax : +91-11-2553 4163, E-mail : info@genaxy.com