


PCR PLATES

AMPLIFY YOUR TRUST
— *NOT JUST* —
NUCLEIC ACID


Genaxy PCR plates are often used by researchers in high-throughput screening thermocycler applications due to their compatibility with most thermal cyclers. The uniform and thin walls enable a precise thermal transfer. The plates can be sealed by the sealing mats, films or foil as per the application.

Unique Features

- » Manufactured with 99.9% pure virgin polypropylene under ISO Class 8 clean room facility
- » Thin, uniform walls to ensure consistent, precise thermal transfer and excellent results
- » Standard and low profile plates are compatible with most major brands of thermal cyclers, like Applied Biosystems[®], BIORAD[®], and Eppendorf[®], Agilent etc.
- » White PCR plates for better results in the analysis of fluorescence signals and consequently with the quantification of replicated DNA.
- » Non-skirt, semi skirt, full skirt and elevated skirt formats are available
- » Printed alphanumeric grid for ease of well identification
- » Certified DNase/RNase, PCR inhibitor, Human DNA and Pyrogen free
- » USFDA 21 CFR Compliant



Genaxy PCR Plates



LOW PROFILE (0.1ML) PLATE



HIGH PROFILE (0.2ML) PLATE



NON-SKIRTED PLATE



NON-SKIRTED ELEVATED WELL PLATE



SEMI-SKIRTED PLATE



ELEVATED SKIRT PLATE



FULL-SKIRTED PLATE



FOR MORE PRODUCT
INFORMATION
SCAN THE QR CODE



PCR Plates

Cat. Number	Description	Qty. / Pack
GEN-96-010-PCR-NS-C	0.1ml, Clear, 96 Well PCR Plate, No Skirt, Low Profile	50
GEN-96-010-PCR-NS-W	0.1ml, White, 96 Well PCR Plate, No Skirt, Low Profile	50
GEN-96-020-PCR-NS-C	0.2ml, Clear, 96 Well PCR Plate, No Skirt, High Profile	50
GEN-96-010-PCR-FS-C	Dual Molded, 0.1ml, Clear, 96 Well PCR Plate, Full Skirt (BioRad Compatible Plate)	50
GEN-96-020-PCR-NSE-C	0.2ml, Clear, 96 Well PCR Plate, No Skirt, elevated wells, High Profile	50
GEN-96-020-PCR-SS-C	0.2ml, Clear, 96 Well PCR Plate, Semi Skirt, High Profile	50
GEN-96-020-ES-ABC	0.2ml, Frosted, 96 Well PCR Plate, Elevated Skirt, High Profile	50
GEN-96-010-ES-ABC	0.1ml, 96-well PCR Plate, fit ABI, Natural, Low Profile	50
GEN-96-010-HS-W-RC	0.1ml, 96-well PCR Plate, PP, Half Skirt, Fit Roche Light Cyler, White, No Sealing Film, Low Profile	50
GEN-96-010-HS-W-SEAL-RC	0.1ml, 96-well PCR Plate PP, Half Skirt, Fit Roche Light Cyler, White, With Sealing Film, Low Profile	50
GEN-384-PCR-FS-C	384 well PCR Plate, Full Skirted, Natural	50
GEN-384-004-FS-W-RC	40µl, 384-well PCR Plate, White Frame, A24/P24 Cut, Fit Roche Light Cyler, White Tubes	50
GEN-384-004-FS-C-ABI	40µl, 384-well PCR Plate, ABI, White Frame, A24 Cut, Clear Tubes	50

Genaxy Scientific Pvt. Ltd.