

# GENAXY HORIZONTAL SYSTEM

## MINI Horizontal Unit

With Genaxy Scientific GENECO range of gel electrophoresis units, you can perform analytical & preparative separations of nucleic acids simply & efficiently. Each system has a wide variety of comb types & comes with a typically waterproof gel casting system. To fulfill your needs, we've created electrophoresis cells that are remarkably advanced, versatile & user-friendly.



### GENECO-H-MINI

## Technical Specification

Model	GENECO-H-MINI
Principal Material	Acrylic
Dimensions (L X B X H)	20 X 10 X 8 cm
No. of Samples	8
Gel Size	10 x 7 cm
Trays (L x W)	10 X 7cm (1), 5 X 7cm (2)
Combs	8 Well x 1.5mm (2) – Acrylic 3 Well Preparative x 3mm- Acrylic Combs
Gel Casting Tray	1 Universal
Platinum Electrode	Red & Black, 1 each
Connecting Cord	Red & Black, 1 each
Lid	1 No.

## Horizontal System Includes

Basic Unit (1 No.)

Platinum Electrode Assembly Removable (2 Nos.)

Lid (1 No.)

8 well comb 1.5 mm (2Nos.)

Gel Running Tray 10 X 7 cm ( 1 No.)

Gel Running Tray 5 X 7 cm ( 2 Nos.)

Connecting Cord (2 Nos.)

3 Well Preparative Comb 3.0 mm

### Accessories

S No.	Catalog No.	Description
1	GENECO-H-MINI-5	Gel Running tray 5 x 7 cm ( 2/pack)
2	GENECO-H-MINI-10	Gel Running tray 10 x 7 cm ( 1/pack)
3	GENECO-H-MINI-GS10	Gel Scoop for 10 X 7 cm gel
4	GENECO-H-LT-15	Levelling table for gel casting 15 x 15 cm
5	GENECO-H-MINI-EL	Spare electrode 2 Nos.
6	GENECO-H-MINI-CC	Connecting Cord

### Combs

S No.	Code	Description
1	GENECO-H-MINI-1.5-8	Acrylic Comb MINI 8 Well 1.5mm
2	GENECO-H-MINI-3.0-3	Acrylic Comb MINI 3 Well 3 mm
3	GENECO-H-MINI-CS	Comb Stand

### Genaxy Scientific Pvt. Ltd

155-156, C-Block, Community Centre, Plot 32 to 57, Behind Janak Cinema, Janak Puri, New Delhi-110 058

Phones : +91-11-4761 9999 (100 Lines) Fax : +91-11-2553 4163

E-mail : info@genaxy.com