



PCR TUBES

AMPLIFY YOUR TRUST
—— *NOT JUST* ——
NUCLEIC ACID

Genaxy PCR tubes are small tubes made of high-quality virgin polypropylene with a conical bottom and unique cap. The most common application of PCR tubes includes various molecular biology biochemistry experiments requiring Polymerase Chain Reaction (PCR) and real-time PCR. As these tubes have ultra-thin walls, they are perfect for accurate and steady thermal transfer in an array of thermocyclers.

Unique Features

PCR SINGLE TUBES


- » Made of medical grade virgin polypropylene conforms to USFDA 21 CFR
- » Manufactured under ISO 8 clean room facility
- » Ultra-smooth, uniform, thin-walled tubes provide better heat transfer
- » Unique cap design with flat cap and domed cap options
- » Clear caps for Real time PCR and frosted cap for labelling are available in 0.2ml tubes
- » Extremely smooth tube bottoms for 100% visibility
- » Certified DNase/RNase, PCR inhibitor, Human DNA, Pyrogen and Phthalate free

PCR STRIP TUBES

- » Medical grade virgin polypropylene 8 strip PCR tubes and caps conforms to USFDA 21 CFR
- » Manufactured under ISO 8 clean room facility
- » Ultra-smooth, uniform, thin walled polypropylene tube-strips for optimum PCR and real time results
- » Unique strip cap design with flat cap and domed cap options (in 0.2ml PCR strip) to prevent sample evaporation
- » Clear well and white well tubes are available for general PCR and qPCR reactions respectively in 0.1ml tube strips
- » 0.1ml 4-tube & 4-cap strips are specifically designed for Rotor-Gene
- » Certified DNase/RNase, PCR inhibitor, Human DNA and Pyrogen free



Genaxy PCR Tubes



LOW PROFILE (0.1ML) STRIP TUBES



HIGH PROFILE (0.2ML) STRIP TUBES



FOR MORE PRODUCT
INFORMATION
SCAN THE QR CODE



PCR Single Tubes

Cat. Number	Description	Qty. / Pack
GEN-01-C-PCR	PCR Tubes, Flat Cap, Non frosted, Clear, 0.1 ml	1000
GEN-02-C-PCR	PCR Tubes, Flat Cap, Non frosted, Clear, 0.2 ml	1000
GEN-02-C-PCR-G	PCR Tubes, Flat Cap, Frosted, Clear, 0.2 ml	1000
GEN-02-DC-PCR	PCR Tubes, Domed Cap, Clear, 0.2 ml	1000
GEN-05-C-PCR	PCR Tubes, Flat Cap, Clear, 0.5 ml	1000

PCR Strip Tubes

Cat. Number	Description	Qty. / Pack
GEN-PT1-1000	0.1ml, 4-Tube & 4-Cap Strips, Roto-Gene [®] Compatible.	250 Strips
GEN-0108-LP-RT-C	8-Strip PCR Tubes with Cap, 0.1 ml, Low Profile, Clear.	120 Strips
GEN-0108-LP-RT-W	8-Strip PCR Tubes with Cap, 0.1 ml, Low Profile, White Tube	120 Strips
GEN-0208-DCC-PCR	8-Strip PCR Tubes with Doomed caps, Clear, 0.2 ml	125 Strips
GEN-0208-FCC-PCR	8-Strip PCR Tubes and Flat Caps, 0.2 ml, Suitable for Optical qPCR	125 Strips
GEN-0208-T-PCR	8-Strip PCR Tubes w/o cap, 0.2 ml	125 Strips
GEN-0208-CPC-PCR	8-Strip, 0.2ml, Domed PCR Caps	125 Strips
GEN-0208-FCP-PCR	8-Strip, 0.2ml, Flat PCR Caps	125 Strips
GEN-0208-RT-PCR	8-Strip, 0.2ml, Optically Clear, Flat Caps for qPCR	125 Strips
GEN-0108-FCCA-PCR	0.1ml, 8-Strip PCR Tubes w/attached Flat caps	120 Strips
GEN-0208-FCCA-PCR	0.2ml, 8-Strip PCR Tubes w/attached Flat caps	120 Strips

Genaxy Scientific Pvt. Ltd.

155-156, C-Block, Community Centre, Plot 32 to 57, Behind Janak Cinema, Janak Puri, New Delhi-110 058
Phones : +91-11-4761 9999 (100 Lines) Fax : +91-11-2553 4163, E-mail : info@genaxy.com